



## GScriptULTRA

### First-Strand Synthesis Kit (Oligo-dT)

Cat. No. MB310-T100

Size: 100 Reactions

Cat. No. MB310-T010

Size: 10 Reactions

Store at -20°C

Minimum shelf life: 1 year

#### Description

The GScriptULTRA First-Strand Synthesis Kit (Oligo-dT) provides a convenient and sensitive of cDNA synthesis from RNA molecules by reverse transcription (RT). The kit consists of five major components: GScriptULTRA RTase, 5X 1st strand buffer, 0.1 M DTT, 50 µM Oligo(dT)<sub>20</sub>, and Random Hexamers primer, and 10 mM dNTP (Deoxyribonucleotide triphosphates) mix. The GScriptULTRA RTase is a recombinant Moloney Murine Leukemia Virus (M-MLV) transcription polymerase expressed in *E. coli* and purified to homogeneity. It has lower RNase H activity and high thermal stability. The enzyme is widely used to synthesize first-strand cDNA at temperatures up to 55°C with increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases. It can generate cDNA from 100 base pairs (bps) to 12 Kilo base pairs (kb). The 5X Reaction Mix buffer is optimized for reverse transcription. The DTT breaks the disulfide bonds, loosen the secondary structure of RNA, and helps in initiation for cDNA synthesis. The Oligo(dT)<sub>20</sub> consist of a stretch of 20 deoxythymidines that anneal to poly(A) tails of eukaryotic mRNAs. The primers are optimized for constructing cDNA libraries. The dNTP solution consists all four nucleotides (dATP, dCTP, dGTP, dTTP), suitable for use in cDNA synthesis.

#### Feature

➤ Time efficiency – included most reagents for reverse transcription reaction.

#### Application

➤ Downstream application for PCR to detection of expressed genes, examination of transcript variants, or generation of cDNA templates for cloning and sequencing.

#### Kit contents

| Contents                        | MB310-T100 | MB310-T010 |
|---------------------------------|------------|------------|
| Oligo(dT) <sub>20</sub> (50 µM) | 100 µl     | 10 µl      |
| 5X 1st strand buffer            | 500 µl     | 50 µl      |
| DTT (0.1 M)                     | 100 µl     | 10 µl      |
| dNTP mix (10 mM)                | 100 µl     | 10 µl      |
| GScriptULTRA RTase              | 40 µl      | 4 µl       |

#### Quality Control

The quality of the GScriptULTRA First-Strand Synthesis Kit (Oligo-dT) is tested on a lot-to-lot basis to ensure consistent product quality.

#### Required Materials

- Microcentrifuge tubes
- RNase-free H<sub>2</sub>O
- PCR instrument or water / Dry bath
- RNase Inhibitor (10 U/µl)

#### GScriptULTRA First-Strand Synthesis Kit (Oligo-dT) Protocol

1. Thaw the GScriptULTRA First-Strand Synthesis Kit (Oligo-dT), template RNA, primers and RNase-free H<sub>2</sub>O on ice. Mix each solution well.
2. Set up the following reaction mixture, and reaction cocktails can be made when multiple reactions are being assembled.

| Components   | Volume (µl) |
|--|-------------|
| Template RNA (10 ng~5 µg total RNA or 10 pg~500 ng mRNA) | Variable    |
| Oligo(dT) <sub>20</sub> (50 µM)                          | 1 µl        |
| 10 mM dNTP Mix   | 1 µl        |
| Add RNase-free H <sub>2</sub> O to                       | 13.6 µl     |

3. Heat for 3-5 minutes at 65°C. Spin briefly and place promptly on ice.
4. Add the following component to the same tube:

| Components                              | Volume (µl) |
|---|-------------|
| 5X 1st strand buffer                    | 4 µl        |
| DTT (0.1 M)                             | 1 µl        |
| RNase Inhibitor (10 U/µl) (recommended) | 1 µl        |
| GScriptULTRA RTase                      | 0.4 µl      |
| Final volume to                         | 20          |

Note: If generating cDNA longer than 5 kb at temperatures above 50°C using a gene-specific primer (GSP) or oligo(dT)<sub>20</sub>, the amount of GScriptULTRA RTase may be raised to 2 µl to increase yield.

5. Incubate at 50°C for 30-60 minutes. Increase the reaction temperature to 55°C for GSP and difficult templates or templates with high secondary.
6. Inactivate enzyme at 70°C for 15 minutes.
7. Store products at -20°C or proceed to PCR. Using 2 µl first-strand cDNA synthesis reaction mixtures. Amplification of some PCR targets (> 1 kb) may require the removal of RNA complementary to the cDNA. To remove RNA complementary to the cDNA, add 1 µl (2 units) of *E. coli* RNase H and incubate at 37°C for 20 minute.

## Troubleshooting

Refer to the table below to troubleshoot problems that you may encounter when you did cDNA synthesis with the kit.

| Problem                   | Possible cause                              | Possible solution  |
|---------------------------|---|--|
| Low yield of PCR products | Incomplete concentration of start materials | Use the appropriate method for the RNA preparation based on the amount of the starting materials.                  |
|                           | RNA degraded                                | Avoid repeated freeze / thaw cycles of the sample.   |
|                           |   | Keep DNA preparations on ice or frozen in order to avoid the degradation.  |
|                           | RNase contaminant                           | Clean everything, use barrier tips, wear gloves and a lab coat, and use RNase-free enzymes, e.g.: RNase inhibitor. |

## Caution:

- During operation, always wear a lab coat, disposable gloves, and protective equipment.
- All products are for research use only .

## Related ordering information

| Cat. No.   | Description  | Size          |
|------------|--|---------------|
| MB101-0500 | <i>Taq</i> DNA Polymerase                              | 500 units     |
| MB203-0100 | OnePCR™  | 100 Reactions |
| DM010-R500 | 1 Kb DNA Ladder RTU                                    | 500 µl        |
| NA017-0100 | Total RNA Isolation Kit (Blood/ Cultured Cell/ Fungus) | 100 Reactions |
| NA020-0100 | Total RNA Isolation Kit (Plant)                        | 100 Reactions |
| NA021-0100 | Total RNA Isolation Kit (Tissue)                       | 100 Reactions |
| RI001-2500 | RiboIN RNase Inhibitor                                 | 2500 U        |