

## MBead Virus Nucleic acid Kit

Cat No. SN011-0100      Size: 100 Reactions  
 Cat No. SN011-0004      Size: 4 Reactions  
 Sample: Up to 300 µl of the virus sample  
 Format: Magnetic Bead System  
 Operation time: Within 10-15 minutes  
 Release volume: 200 µl



### Description

This MBead Virus Nucleic acid Kit is designed specifically for the simultaneous virus DNA/RNA purification from the plasma, serum, body fluid or supernatant of virus-infected cell cultures. Its unique buffer system will efficiently lyse cells and degrade proteins, allowing for the nucleic acid to be easily bound by the surface of the magnetic beads. The other non-specific binding particles are removed with a wash buffer, and the nucleic acid is released into the Release Buffer. The nucleic acid can be purified manually within 10-15 minutes (using most magnetic separators) or the kit can be easily adapted to satisfy most automated nucleic acid purification systems.

### Features

- Fast, reproducible, and easy processing by using a magnetic bead system.
- To isolate high quality nucleic acid.
- Isolated nucleic acid is compatible with various downstream applications.

### Applications

- Restriction enzyme digestion.
- PCR amplification.
- Real-Time PCR assay.

### Kit Contents

Contents	SN011-0100	SN011-0004
Magnetic Bead	2 ml X 1 vial	80 µl X 1 vial
Lysis Buffer	30 ml X 1 bottle	1.5 ml X 1 vial
Wash Buffer	80 ml X 1 bottle	2 ml X 2 vials
Release Buffer	20 ml X 1 bottle	1 ml X 1 vial

### Quality Control

The quality of the MBead Virus Genomic DNA Kit is tested on a lot-to-lot basis to ensure consistent product quality.

### Required Materials

- Absolute ethanol
- 1.5 ml microcentrifuge tubes
- Magnetic separator
- Water bath/ Dry bath

### MBead Virus Nucleic acid Kit Protocol

#### Step 1 Lysis

1. Transfer the virus sample (up to 300 µl) into a 1.5 ml microcentrifuge tube and add 300 µl of the Lysis Buffer.
2. Mix well and incubate at 65°C for 5 minutes. During this time, pre-heat the Release Buffer to 65°C for the Step 4.
3. Add 300 µl of the absolute ethanol to the lysate and mix well.

#### Step 2 DNA Binding

1. Add 20 µl of the Magnetic Beads. Mix well by gently shaking for 3 minutes.
2. Place the tube in a magnetic separator for 30 seconds.
3. Remove the solution (If the mixture becomes viscous, increase magnetic bead separation time).


#### Step 3 Wash

1. Add 800 µl of the Wash Buffer and mix well (Following the wash, the mixture will no longer be viscous).
2. Place the tube in a magnetic separator for 30 seconds. Remove the solution.

#### Step 4 Release

1. Add 200 µl of the Release Buffer (pre-heated to 65°C) and mix well.
2. Incubate for 3 minutes at 65°C (During the incubation, shake the tube vigorously every minute).
3. Place the tube in a magnetic separator for 1 minute.
4. Carefully transfer ONLY the clean portion of the solution to a clean tube.

**NOTE:** Be sure and allow the magnetic beads to disperse completely during the binding, wash and release steps.



### Troubleshooting

Refer to the table below to troubleshoot problems that you may encounter when you did DNA/ RNA isolation with the kit.

Problem	Cause	Solution
Low yields of nucleic acid	Low virus recovery	Optimize virus production.
	Release rate was too fast	Following the release step 4, incubate for 3 minutes at 65°C.
High amount of cellular DNA	Excess of cellular DNA	Decrease the input amount of viral supernatant.

## Related Ordering Information

Cat. No.	Description	Size
SM101-0500	<i>Taq</i> DNA polymerase	500 U
SM200-0100	PCR SUPERMIX	100 Reactions
SM201-0100	Hot Start SUPERMIX	100 Reactions
SA001-0500	AGAROSE Tablet, 0.5g	100 Tablets
SL001-1000	Novel Juice Supplied in 6X Loading Buffer	1 ml
SD003-R600	100 bp DNA Ladder H3 RTU	600 µl
SD010-R600	1 Kb DNA Ladder RTU	600 µl
SD013-R600	XLarge DNA Ladder RTU	600 µl
ST040-4000	100 mM dNTP Set	4 x 1 ml
ST046-1000	100 mM dNTP Set	4 x 250 µl
ST025-1000	2.5 mM dNTP Mix	1 ml
ST010-1000	10 mM dNTP Mix	1 ml

## Caution

- Check buffers before use for salt precipitation. Re-dissolve any precipitate by warming up to 37°C.
- During operation, always wear a lab coat, disposable gloves, and protective equipment.
- All products are for research use only.