

Agrobacterium (AB) Minimal Medium

Agrobacterium (AB) minimal medium contains AB salts, AB buffers, and glucose. It is a minimal growth medium that can be used in bacterial cultures. Compared to M9 medium.

Components	g/L
Ammonium Chloride	1.0
Magnesium Sulfate, Heptahydrate	0.3
Potassium Chloride	0.15
Calcium Chloride	0.01
Ferrous Sulfate, Heptahydrate	0.00075
Potassium Phosphate, Dibasic	3.0
Sodium Dihydrogen Phosphate	1.0
MES Free Acid	0.0195
Glucose	5.0

Storage

Store at room temperature. Store prepared medium at 2-8°C.