

Chu (N6) Medium w/ Vitamins

Working concentration: 4.1 g/L

Micronutrients	mg/L
Boric acid	1.6
Manganese Sulfate•H ₂ O	3.3
Potassium Iodide	0.8
Zinc Sulfate•7H ₂ O	1.5
EDTA Disodium Salt Dihydrate	37.25
Ferric Sulfate•7H ₂ O	27.85
Macronutrients	
Calcium Chloride, anhydrous	125.33
Magnesium Sulfate, anhydrous	90.7
Ammonium Sulfate	463
Potassium Nitrate	2830
Potassium Phosphate Monobasic	400
Vitamins	
Thiamine HCl	1
Pyridoxine HCl	0.5
Nicotinic Acid	0.5
Glycine	2
Total	4100

