

## リムーバブルニードル (RN) 交換針

| 針先<br>ゲージ | 2.5-100 $\mu$ L シリンジ用 |             |              |
|-----------|-----------------------|-------------|--------------|
|           | PT2<br>(鋭角)           | PT3<br>(直角) | PT5<br>(横穴型) |
| 26s       | 7758-02               | 7768-01     | 7784-07      |
| 26        | 7758-04               | 7768-02     | 7784-08      |
| 22s       | 7758-03               | 7770-01     | 7784-05      |
| 22        | 7758-01               | 7770-02     | 7784-06      |
| 25s       | —                     | 8647-01     | —            |

| 250 $\mu$ L-100mL シリンジ用 |             |              |
|-------------------------|-------------|--------------|
| PT2<br>(鋭角)             | PT3<br>(直角) | PT5<br>(横穴型) |
| 7779-02                 | 7780-01     | 7784-03      |
| 7779-04                 | 7780-02     | 7784-04      |
| 7779-03                 | 7780-03     | 7784-01      |
| 7779-01                 | 7780-04     | 7784-02      |
| —                       | 8648-01     | —            |

## Kel-F ハブ付き交換針 (KFタイプ)

| ゲージ | 針先     |        | Model  | PT2<br>(鋭角) | PT3<br>(直角) | PT5<br>(横穴型) |
|-----|--------|--------|--------|-------------|-------------|--------------|
|     | 外径(mm) | 内径(mm) |        |             |             |              |
| 33  | 0.21   | 0.11   | —      | —           | —           | —            |
| 32  | 0.24   | 0.11   | —      | —           | —           | —            |
| 31  | 0.26   | 0.13   | KF731  | 90131       | 90531       | —            |
| 30  | 0.31   | 0.16   | KF730  | 90130       | 90530       | —            |
| 29  | 0.34   | 0.18   | KF729  | 90129       | 90529       | —            |
| 28  | 0.36   | 0.18   | KF728  | 90128       | 90528       | —            |
| 27  | 0.41   | 0.21   | KF727  | 90127       | 90532       | —            |
| 26s | 0.47   | 0.13   | KF726s | 90139       | 90539       | 7746-12      |
| 26  | 0.46   | 0.26   | KF726  | 90126       | 90533       | 7746-10      |
| 25  | 0.51   | 0.26   | KF725  | 90125       | 90525       | 7746-09      |
| 24  | 0.57   | 0.31   | KF724  | 90124       | 90524       | 7746-08      |
| 23  | 0.64   | 0.34   | KF723  | 90123       | 90523       | 7746-07      |
| 22s | 0.72   | 0.15   | KF722s | 90138       | 90534       | 7746-11      |
| 22  | 0.72   | 0.41   | KF722  | 90122       | 90134       | 7746-06      |
| 21  | 0.82   | 0.51   | KF721  | 90121       | 90521       | 7746-05      |
| 20  | 0.91   | 0.60   | KF720  | 90120       | 90520       | 7746-04      |
| 19  | 1.07   | 0.69   | KF719  | 90119       | 90519       | 7746-13      |
| 18  | 1.27   | 0.84   | KF718  | 90118       | 90535       | 7746-03      |
| 17  | 1.47   | 1.07   | KF717  | 90117       | 90517       | 7746-02      |
| 16  | 1.65   | 1.19   | KF716  | 90116       | 90516       | 7746-01      |
| 15  | 1.83   | 1.37   | KF715  | 90115       | 90515       | 7746-14      |

## メタルハブ付き交換針 (Nタイプ)

| Model | 針先          |             |              |
|-------|-------------|-------------|--------------|
|       | PT2<br>(鋭角) | PT3<br>(直角) | PT5<br>(横穴型) |
| N733  | 90033       | 91033       | —            |
| N732  | 90032       | 91032       | —            |
| N731  | 90031       | 91031       | —            |
| N730  | 90030       | 91030       | —            |
| N729  | 90029       | 91029       | —            |
| N728  | 90028       | 91028       | —            |
| N727  | 90027       | 91027       | —            |
| N726s | 90039       | 91039       | 7729-01      |
| N726  | 90026       | 91026       | 7729-03      |
| N725  | 90025       | 91025       | 7729-04      |
| N724  | 90024       | 91024       | 7729-05      |
| N723  | 90023       | 91023       | 7729-06      |
| N722s | 90038       | 91038       | 7729-02      |
| N722  | 90022       | 91022       | 7729-07      |
| N721  | 90021       | 91021       | 7729-08      |
| N720  | 90020       | 91020       | 7729-09      |
| N719  | 90019       | 91019       | 7729-14      |
| N718  | 90018       | 91018       | 7729-10      |
| N717  | 90017       | 91017       | 7729-11      |
| N716  | 90016       | 91016       | 7729-12      |
| N715  | 90015       | 91015       | 7729-13      |