

Create Date: 2017-03-28 00:00:00/ Revision Date: 2019-05-07 01:00:00/ Print Date: 2019-10-29

---

**SECTION 1: COMPANY AND PRODUCT INFORMATION**

---

**1.1 Product identifiers**

Product name : Crystal Violet 0.01%  
Product code : 40140028  
CAS number : Not Available  
Synonyms : Not Available

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Identified uses : For research and laboratory use only.

**1.3 Details of the supplier of the safety data sheet**

Address : Genelix International Inc, dba bioWORLD  
4150 Tuller Rd. Suite 228  
Dublin, OH 43017

Email : tech@bio-world.com

Phone : 614-792-8680, Toll free: 1-888-bio-PLUS

Fax : 614-792-8685

**1.4 Emergency telephone number**

Emergency phone : 1-888-bio-PLUS

---

**SECTION 2: HAZARDS IDENTIFICATION**

---

**2.1 Classification of substance or mixture**

Eye irritation (Category 2A), H319  
Carcinogenicity (Category 2), H351  
Acute aquatic toxicity (Category 2), H401  
Chronic aquatic toxicity (Category 2), H411

**2.2 GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Eye irritation (Category 2A), H319  
Carcinogenicity (Category 2), H351  
Acute aquatic toxicity (Category 2), H401  
Chronic aquatic toxicity (Category 2), H411

**2.3 Label elements and precautionary statements**

Pictogram : 

Signal word : Warning

Hazard statement(s) : H319 - Causes serious eye irritation.  
H351 - Suspected of causing cancer.  
H411 - Toxic to aquatic life with long lasting effects.

Precautionary statement(s) : P201 - Obtain special instructions before use.  
P202 - Do not handle until all safety precautions have been read and understood.  
P264 - Wash skin thoroughly after handling.  
P273 - Avoid release to the environment.  
P280 - Wear protective gloves/eye protection/face protection.  
P305+351+338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P308+313 - IF exposed or concerned: Get medical advice/attention.  
P337+313 - If eye irritation persists: Get medical advice/attention.  
P391 - Collect spillage.  
P405 - Store locked up.  
P501 - Dispose of contents/container to an approved waste disposal plant.

## 2.4 Hazards not otherwise classified (HNOC) or not covered by GHS

No unclassified hazards known.

## 2.5 NFPA Rating

Health hazard : 2

Fire hazard : 0

Reactivity hazard : 0

## 2.6 HMIS Rating

Health hazard : 2

Chronic health hazard : \*

Reactivity hazard : 0

Flammability : 0

Physical hazard : 0

---

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

---

### 3.1 Substances

Substance	CAS#	EC#	Concentration
Crystal Violet M.F: C <sub>25</sub> H <sub>30</sub> N <sub>3</sub> Cl M.W: 407.98 g/mol	548-62-9	208-953-6	<= 1%

### 3.2 Hazardous components & classification

Eye Irrit. 2A; Carc. 2; Acute Aqua. Tox. 2; Chronic Aqua. Tox. 2; H319, H351, H401, H411.

---

## SECTION 4: FIRST AID MEASURES

---

### 4.1 Description of first aid measures

#### General advice

Consult a physician if symptoms are severe or persistent. Provide this data sheet to medical personnel. If product is spilled or leaked, evacuate area.

#### In case of inhalation

If inhaled, move person to fresh air and monitor breathing. If not breathing, give artificial ventilation. Consult a physician if symptoms are severe or persistent.

#### **In case of skin contact**

Immediately wash with excess soap and water. If spilled on clothing, remove all affected clothing. Consult a physician if symptoms are severe or persistent.

#### **In case of eye contact**

Flush eyes with water or eye wash solution as a precaution for 15 minutes. Consult a physician if symptoms are severe or persistent.

#### **In case of ingestion**

Only induce vomiting if recommended by medical personnel. If subject is unconscious, do not give anything by mouth. If conscious, rinse mouth with water. Consult a physician if symptoms are severe or persistent.

#### **4.2 Most important symptoms and effects, both acute and delayed**

All known important symptoms are described in Section 2 and/or Section 11. No other important symptoms to report.

#### **4.3 Indication of any immediate medical attention and special treatment needed**

No special treatment indicated. Provide treatment in accordance with exhibited systems.

---

### **SECTION 5: FIREFIGHTING MEASURES**

---

#### **5.1 Suitable extinguishing media**

Water spray, alcohol-resistant foam, dry chemical, and carbon dioxide extinguishers are suitable.

#### **5.2 Unsuitable extinguishing media**

No known unsuitable extinguishing media.

#### **5.3 Special hazards arising from the substance**

No special hazards are known to exist regarding combustion of this product.

#### **5.4 Advice for firefighters**

Wear protective gear, such as self-contained breathing apparatus, if necessary.

---

### **SECTION 6: ACCIDENTAL RELEASE MEASURES**

---

#### **6.1 Personal precautions, protective equipment and emergency procedures**

Provide suitable ventilation. Use any necessary personal protective equipment. Avoid contact with skin and eyes, and avoid creation and inhalation of vapor or dust. Keep all unnecessary personnel away.

**For personal protection see section 8**

#### **6.2 Environmental precautions**

Prevent product from entering public sewers and waterways.

#### **6.3 Methods and material for containment and cleaning up**

Use inert absorbent material to absorb any spilled or leaked product and dispose of as hazardous waste. Keep in suitable, closed containers for disposal.

**For proper disposal see section 13**

---

### **SECTION 7: HANDLING AND STORAGE**

---

#### **7.1 Precautions for safe handling**

Provide suitable ventilation. Wear any necessary personal protective equipment. Avoid inhalation of vapor or mist.

**For precautions see section 2**

## 7.2 Conditions for safe storage, including any incompatibilities

Storage conditions : Store upright, closed container in arid, ventilated environment.

Incompatible materials : Not Available

---

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

---

### 8.1 Control parameters

This product is not known to contain any substances with occupational exposure limit values.

### 8.2 Engineering controls

Follow good industrial hygiene and safety practices when handling product.

### 8.3 Personal protective equipment

Eye/face protection : Use only government-approved safety glasses with side-shields.

Skin protection : Use gloves when handling product. Inspect gloves before use to ensure suitability for use. Remove without exposing skin to the gloves outer surface. Discard used gloves according to all pertinent laws and/or current good practices (cGXP). Wash hands with soap and water.

Body protection : Wear appropriate clothing. Ensure clothing is in good condition, with no holes or tears. When selecting clothing, consider the concentration and amount of substance to be handled.

Respiratory protection : Use only approved respirators and components which comply with CDC and NIOSH (US) or CEN (EU) regulations. Required only when vapors or aerosols are created.

Control of environmental exposure : Prevent product from entering the environment, especially through public sewers or waterways.

General hygiene considerations : Comply with general industrial hygiene practice guidelines.

---

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

---

### 9.1 Information on basic physical and chemical properties

- a) Appearance : Dark Violet
- b) Physical states : Liquid
- c) Odor : Not available
- d) Odor threshold : Not available
- e) Melting point : Not available
- f) Boiling point range : Not available
- g) pH : Not available

h)	Density	: Not available
i)	Flash point	: Not available
j)	Evaporation rate	: Not available
k)	Flammability	: Not available
l)	Upper/lower flammability or explosive limits:	: Not available
m)	Vapor pressure	: Not available
n)	Vapor density	: Not available
o)	Relative density	: Not available
p)	Water solubility	: Not available
q)	Partition coefficient:n-octanol/water	: Not available
r)	Autoignition temperature	: Not available
s)	Decomposition temperature	: Not available
t)	Kinematic viscosity	: Not available
u)	Explosive properties	: Not available
v)	Oxidizing properties	: Not available
w)	Solubility in other solvents	: Not available
x)	Surface tension	: Not available

---

## SECTION 10: STABILITY AND REACTIVITY

---

### 10.1 Reactivity

No special reactivity is known.

### 10.2 Chemical stability

Product is stable when stored and used as recommended.

### 10.3 Stability note(s)

No special or unusual instability known.

### 10.4 Polymerization

No known polymerization.

### 10.5 Possibility of hazardous reactions

No hazardous reactions are known.

### 10.6 Incompatible materials

Not available.

### 10.7 Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

---

## SECTION 11: TOXICOLOGICAL INFORMATION

---

### 11.1 Acute toxicity

LD50 Oral : No toxicity data available.

LD50 Dermal : No toxicity data available.

LC50 Inhalation : No toxicity data available.

**11.2 Skin corrosion/irritation**

No skin/corrosion irritation data available.

**11.3 Serious eye damage/eye irritation**

No eye damage/irritation data available.

**11.4 Respiratory or skin sensitization**

No sensitization data available.

**11.5 Germ cell mutagenicity**

No mutagenicity data available.

**11.6 Carcinogenicity**

IARC : Product and components are not regulated by the IARC.

ACGIH : Product and components are not regulated by the ACGIH.

NTP : Product and components are not regulated by the NTP.

OSHA : Product and components are not regulated by OSHA.

**11.7 Reproductive toxicity**

No reproductive toxicity data available.

**11.8 Specific target organ toxicity – single exposure**

No specific organ toxicity data available.

**11.9 Specific target organ toxicity – repeated exposure**

No specific organ toxicity data available.

**11.10 Aspiration hazard**

No aspiration hazard data available.

**11.11 Additional Information**

RTECS: No RTECS data available.

---

**SECTION 12: ECOLOGICAL INFORMATION**

---

**12.1 Toxicity**

No ecological toxicity data available.

**12.2 Aquatic toxicity**

No aquatic toxicity data available.

**12.3 Persistence and degradability**

No persistence/degradability data available.

**12.4 Bioaccumulative potential**

No bioaccumulation data available.

**12.5 Mobility in soil**

No soil mobility data available.

**12.6 Results of PBT and vPvB assessment**

PBT/vPvB assessment not available as chemical safety assessment is not required/was not conducted.

## 12.7 Other adverse effect

Toxic to aquatic life with long lasting effects.

---

## SECTION 13: DISPOSAL CONSIDERATIONS

---

### 13.1 Product

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. Chemical waste generators must also consult and adhere to local, regional, and national hazardous waste regulations to ensure complete and accurate classification. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

### 13.2 Packaging

Packaging should be disposed of in the same manner as unused product.

### 13.3 Recommendation

Disposal must be made according to official regulations.

---

## SECTION 14: TRANSPORTATION INFORMATION

---

### 14.1 DOT (US)

Not a dangerous good under DOT(US) regulations.

### 14.2 IMDG

UN# 3082

Class: 9

Packing Group: III

EMS No: F-A, S-F

Proper Shipping Name: Environmentally hazardous substance, liquid, n.o.s. (Crystal Violet)

### 14.3 IATA

UN# 3082

Class: 9

Packing Group: III

EMS No: F-A, S-F

Proper Shipping Name: Environmentally hazardous substance, liquid, n.o.s. (Crystal Violet)

EHS-Mark required (ADR 2.2.9.1.10, IMDG code 2.10.3) for single packagings and combination packagings containing inner packagings with Dangerous Goods > 5L for liquids or > 5kg for solids.

---

## SECTION 15: REGULATORY INFORMATION

---

### 15.1 SARA

SARA 302: This product and components are not subject to the reporting requirements of SARA Title III, Section 302.

SARA 313: This product does not contain any components with known CAS numbers that exceed the threshold reporting levels established by SARA Title III, Section 313.

SARA 311/312: Acute Health Hazard, Chronic Health Hazard

### 15.2 Clean water act (CWA)

No chemicals are present in this product that are subject to regulation under the Clean Water Act.

### 15.3 Right to know components

Massachusetts	: No chemicals are present which require disclosure under the Massachusetts Right to Know Act.
Pennsylvania	: No chemicals are present which require disclosure under the Pennsylvania Right to Know Act.
New Jersey	: C.I. Basic Violet 3 CAS No.: 548-62-9
California proposition 65 components	: This product contains no chemicals which are known to the State of California to cause cancer, or birth defects or other reproductive harm.

---

## **SECTION 16: OTHER INFORMATION**

---

### **16.1 Disclaimer**

This product is offered by bioWORLD for research, laboratory or further manufacturing use. Not for human use or consumption. The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchant-ability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall bio-WORLD.com be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if bio-WORLD.com has been advised of the possibility of such damages.

### **16.2 Copyright**

Copyright 2019 www.bio-WORLD.com. License granted to make paper copies for internal use only.

### **16.2 Preparation Information**

Create Date: 2017-03-28 00:00:00/ Revision Date: 2019-05-07 01:00:00/ Print Date: 2019-10-29