

Gamborg B5 Medium

Catalog Number:

SKU#

30630028

Description: Contains micronutrients and macronutrients as described by Gamborg, et al. (1968).

Media
Concentration: 3.1 g/L

Components (mg/L)	
Ammonium Sulfate	134
Boric Acid	3
Calcium Chloride, Anhydrous	113.24
Cobalt Chloride, Hexahydrate	0.025
Cupric Sulfate, Pentahydrate	0.025
Disodium EDTA, Dihydrate	37.26
Ferrous Sulfate, Heptahydrate	27.8
Magnesium Sulfate, Anhydrous	122.09
Manganese Sulfate, Monohydrate	10
Molybdic Acid Sodium Salt, Dihydrate	0.25
Potassium Iodide	0.75
Potassium Nitrate	2500
Sodium Phosphate Monobasic	150
Zinc Sulfate, Heptahydrate	2



plantmedia
a division of bioWorld
www.plantmedia.com
1.888.bio.PLUS

USA

fax: 614-792-8685

4150 Tuller Road Suite 228 Dublin, Ohio 43017

