



Health	2
Fire	0
Reactivity	0
Personal Protection	

Material Safety Data Sheet

TBS buffer 10X

1: Chemical Product and Company Identification

Product Name: TBS buffer 10X

Contact information:

Catalog Code: 42020056-1, 42020056-2, 42020056-3

bioWORLD

4150 Tuller Rd

Suite 228, Dublin, OH 43017, USA

CAS#: Mixture.

For information, call: 614-792-8680

Fax Number: 614-792-8685

Order Online: www.bio-world.com

RTECS: Mixture.

TSCA: TSCA 8(b) inventory : Components were listed.

E-mail : save@bio-world.com

CI#: Not available.

Synonym: Not available.

In case of a chemical emergency, contact local hospital.

Chemical Name: Tris buffered saline

For non-emergency assistance, call: 1-800-860-9729

Chemical Formula: Mixture.

2: Composition and Information on Ingredients

Composition:

Name	CAS#	% by Weight
Tris	77-86-1	<4
Sodium Chloride	7647-14-5	<5
Potassium Chloride	7447-40-7	<1
Water	7732-18-5	>90

Toxicological Data on Ingredients:

Sodium chloride

LC₅₀ (Inhalation, Rat): >42 gm/m³/1H

LD₅₀ (Oral, Mouse): 4 gm/kg

LD₅₀ (Oral, Rat): 3000 mg/kg

Potassium chloride

LD₅₀ (Oral, Rat): 2600 mg/kg

3: Hazards Identification

Potential Acute Health Effects: Irritant.

Irritating to eyes, respiratory system and skin.

Potential Chronic Health Effects:

CARCINOGENIC EFFECTS: Not available.

MUTAGENIC EFFECTS: Not available.

TERATOGENIC EFFECTS: Not available.

DEVELOPMENTAL TOXICITY: Not available.

4: First Aid Measures

Eye Contact: In case of contact, immediately flush eyes with copious amounts of water for at least 15 minutes.

Skin Contact: In case of contact, immediately wash skin with soap and copious amounts of water.

Serious Skin Contact: Wash with a disinfectant soap and cover the contaminated skin with an anti-bacterial cream. Seek medical attention.

Inhalation: If inhaled, remove to fresh air. If not breathing give artificial respiration. If breathing is difficult, give oxygen.

Serious Inhalation: Not available.

Ingestion: If swallowed, wash out mouth with water provided person is conscious. Call a physician.

Serious Ingestion: Not available.

5: Fire and Explosion Data

Flammability of the Product: Noncombustible.

Auto-Ignition Temperature: Not available.

Flash Points: Not available.

Flammable Limits: Not available.

Product of Combustion: Nature of decomposition products not known.

Fire Hazards in Presence of Various Substances: Not available.

Explosion Hazards in Presence of Various Substances:
Not available.

Fire Fighting Media and Instructions:

SMALL FIRE: Dry chemical powder.

LARGE FIRE: Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.

Special Remarks on Fire Hazards: Not available.

Special Remarks on Explosion Hazards: Not available.

6: Accidental Release Measures

Small Spill: Absorb on sand or vermiculite and place in closed containers for disposal. Ventilate area and wash spill site after material pickup is complete.

Large Spill: Use a shovel to put the material into a convenient waste disposal container. Finish cleaning by spreading water on the contaminated surface and allow to evacuate through the sanitary system.

7: Handling and Storage

Precautions: Do not breathe vapor. Avoid contact with eyes, skin, and clothing. Avoid prolonged or repeated exposure.

Storage: Keep tightly closed. Store at 2-8°C.

8: Exposure Controls/Personal Protection

Engineering Controls: Safety shower and eye bath. Mechanical exhaust required.

Personal Protection: Respiratory : Government approved respirator.

Hand: Chemical resistant gloves.

Eye: Safety glasses.

Personal Protection in Case of a Large Spill: Splash goggles. Full suit. Boots. Dust respirator. Gloves.

Exposure Limits: Not available.

9: Physical and Chemical Properties

Physical State and Appearance: Liquid.

Odor: Not available.

Taste: Not available.

Molecular Weight: Mixture.

Color: Colorless.

pH(1% soln/water): 7.25 - 7.55

Boiling Point: > 100 °C

Melting Point: Not available.

Critical Temperature: Not available.

Specific Gravity: >1

Vapour Pressure: Not available.

Vapor Density: >1

Volatility: Not available.

Odor Threshold: Not available.

Water/Oil Dist. Coeff: Not available.

Ionicity (in Water): Not available.

Dispersion Properties: Not available.

Solubility: Soluble in water.

10: Stability and Reactivity Data

Stability: Stable.

Instability Temperature: Not available.

Conditions of Instability: Not available.

Incompatibility with various substances: Strong oxidizing agents, Bases.

Corrosivity: Not available.

Special Remarks on Reactivity: Not available.

Special Remarks on Corrosivity: Not available.

Polymerization: Will not occur.

11: Toxicological Information

Routes of Entry: Skin contact. Eye contact. Ingestion. Inhalation.

Toxicity to Animals:

Sodium chloride
LC₅₀ (Inhalation, Rat): >42 gm/m³/1H
LD₅₀ (Oral, Mouse): 4 gm/kg
LD₅₀ (Oral, Rat): 3000 mg/kg
LD₅₀ (Skin, Rabbit): >10 gm/kg
Potassium chloride
LD₅₀ (Oral, Mouse): 1500 mg/kg
LD₅₀ (Oral, Rat): 2600 mg/kg

Chronic Effects on Humans: Not available.

Other Toxic Effects on Humans: Skin Contact: Causes skin irritation.

Skin Absorption: May be harmful if absorbed through the skin.

Eye Contact: Causes eye irritation.

Inhalation: May be harmful if inhaled. Material is irritating to mucous membranes and upper respiratory tract.

Ingestion: May be harmful if swallowed.

Special Remarks on Toxicity to Animals: Not available.

Special Remarks on Chronic Effects on Humans: Not available.

Special Remarks on other Toxic Effects on Humans: Not available.

12: Ecological Information

Ecotoxicity: Not available.

BOD5 and COD: Not available.

Product of Biodegradation: Not available.

Toxicity of the Products of Biodegradation: Not available.

Special Remarks on the Products of Biodegradation: Not available.

13: Disposal Considerations

Waste Disposal: Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Observe all federal, state, and local environmental regulations.

14: Transport Information

DOT Classification: Not a DOT controlled material.

Identification: Not available.

Special Provisions for Transport: This substance is considered to be non-hazardous for transport.

15: Other Regulatory Information

Federal and State Regulations: TSCA 8(b) inventory : Tris; Sodium chloride; Potassium chloride; Water

Other Regulations: EU ADDITIONAL CLASSIFICATION

Symbol of Danger: Xi

Indication of Danger: Irritant.

R: 36/37/38

Risk Statements: Irritating to eyes, respiratory system and skin.

S: 26-36

Safety Statements: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Other Classifications:

WHMIS(Canada): Not available.

DSCL(EEC): Not available.

HMIS(U.S.A):

Health Hazard: 2

Fire Hazard: 0

Reactivity: 0

Personal Protection:

National Fire Protection Association (U.S.A):

Health: 2

Flammability: 0

Reactivity: 0

Specific hazard:

Protective Equipment: Gloves. Lab coat. Dust respirator. Be sure to use an approved/certified respirator or equivalent. Safety goggles.

16: Other Information

References: Not available.

Other Special Considerations: Not available.

Created: 21/09/2006 12:07:59 PM

Last Updated: 17/12/2008 1:15:49 PM

Disclaimer

This product is offered by bioWORLD for research, laboratory or further manufacturing use. Not for human use or consumption.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall Bio-World.com be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if Bio-World.com has been advised of the possibility of such damages.