

Catalog Number:

40120830

## Hank's Balanced Salt Solution (HBSS) with HEPES, Sterile

**Description:** Hank's Balanced Salt Solution (HBSS) without phenol red and sodium bicarbonate and with the addition of 10mM HEPES, pH 7.4, Sterile.

**Composition:**

Potassium Chloride	5.4mM
Sodium phosphate dibasic heptahydrate	0.3mM
Potassium phosphate monobasic anhydrous	0.4mM
HEPES	10mM
Calcium Chloride	1.3mM
Magnesium Chloride Hexahydrate	0.5mM
Magnesium Sulfate, Heptahydrate	0.6mM
Sodium Chloride	137mM
D-glucose	5.6mM

*Serving Molecular Life Sciences*

[www.bio-WORLD.com](http://www.bio-WORLD.com)

**bioWORLD**

made in **USA**

4150 Tuller Road, Suite 228 Dublin, Ohio 43017 • p: 1.888.BIOPLUS • f: 1.614.792.8685 • e: [save@bio-WORLD.com](mailto:save@bio-WORLD.com)