

Aspiration Pipette

(#94001/94002/94005/94010)

1. Description

SPL Aspiration Pipettes are useful for quick collection and discardment of liquid samples.

2. Features

- ❖ No graduation scales
- ❖ No plugs

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase free	:	Passed
DNA free	:	Passed

4. Product Range

Cat.No.	Material	Working Volume (ml)	Sterile	Packaging
94001	PS	1.00	+	1 / 200 / 800
94002	PS	2.00	+	1 / 150 / 600
94005	PS	5.00	+	1 / 100 / 400
94010	PS	10.00	+	1 / 100 / 400

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs



SPL-TDS-ASP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 19th, 2016


7. Specifications & Technical Drawing

(1) 1 ml (# 94001)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 277.00 mm Diameter: 4.60 mm Volume: 1.00 ml	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 1 - Pieces/Inbox: 200 - Pieces/Case: 800

- Technical Drawing




(2) 2 ml (# 94002)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 277.00 mm Diameter: 6.40 mm Volume: 2.00 ml	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 1 - Pieces/Inbox: 150 - Pieces/Case: 600

- Technical Drawing




(3) 5 ml (# 94005)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 338.00 mm Diameter: 8.30 mm Volume: 5.00 ml	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 1 - Pieces/Inbox: 100 - Pieces/Case: 400

- Technical Drawing




(4) 10 ml (# 94010)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 338.00 mm Diameter: 9.90 mm Volume: 10.00 ml	Crystal-grade polystyrene (PS)	- Pieces/Sleeve: 1 - Pieces/Inbox: 100 - Pieces/Case: 400

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-ASP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© SPL Life Sciences Co., Ltd., Version 2, July 19th, 2016