

Novel Green (10000X)

1. IDENTIFICATION OF THE PREPARATION/SUBSTANCE AND OF THE COMPANY

Preparation Name: Novel Green (10000X) Catalogue number of the product: SL002-0500, SL002-0010 Use of the preparation: For laboratory use. Company identification: GeneDireX, Inc. Website: www.genedirex.com

2. COMPOSITION/ INFORMATION ON INGREDIENTS

General information: The product contains: Solvents Ingredients: CAS#			
Name	CAS No.	EC No.	Weight %
Novel Green Dye			<0.01% (v/v)
Dimethyl sulfoxide (DMSO)	67-68-5	200-664-3	> 90 % (v/v)

3. HAZARDS IDENTIFICATION

Warning! May cause eye, skin and respiratory tract irritation. Combustible liquid and vapor. Human health hazards: Inhalation High concentration: May cause respiratory system irritation. May cause skin irritation. May be absorbed through the skin. Skin contact May cause temporary eye irritation. Exposed individuals may experience eye Eye contact redness, tearing and discomfort. May be uncomfortable if swallowed. Ingestion The liver and kidneys injury may occur. Chronic exposure Target organ(s) Eye, Skin, Respiratory system, Kidney and Liver. Hazards deriving from physical-chemical properties The product is combustible, but not flammable. This product is hazardous according to OSHA 29CFR 1910.1200. OSHA regulatory status Environmental hazards N.A.

4. FIRST-AID MEASURES

Ingestion:

If swallowed, rinse mouth thoroughly with water and give large amounts of milk or water to people not unconscious. Give medical attention if any discomfort continues.

Inhalation:

If inhaled, move injured person into fresh air and keep person calm under observation.

Give medical attention if any discomfort continues.

Skin contact:

Immediately flush skin with water for at least 15 minutes while removing contaminated clothes.

Get medical attention if irritation persists after washing.

Eye contact:

In case of contact, immediately flush with plenty of water for up to 15 minutes.

Remove any contact lenses and open eyes wide apart. Get medical assistance if irritation occurs.



5. FIRE-FIGHTING MEASURES

Suitable extinguishing media:

Water spray, Carbon Dioxide, dry chemical powder or appropriate foam.

Not suitable extinguishing media:

None.

Special fire fighting procedures:

Use standard firefighting procedures. Containers close to fire should be removed or cooled with water. Unusual fire and explosion hazards:

During fire, gases hazardous to health may be formed.

Hazardous combustion products:

Carbon dioxide, carbon monoxide, formaldehyde, sulfur oxides.

Protective equipment:

Wear self-contained breathing apparatus and protective clothing. Follow the general fire precautions indicated in the workplace.

Flammable class:

NFPA rating fire=2

6. ACCIDENTAL RELEASE MEASURES

Personal precautions:

Wear appropriate gloves (see Section 8), safety glasses to clean small releases/spillage.

In case of vapor and spray mist, wear suitable respiratory protection means. Environmental precautions: Avoid discharge into drains, water courses or onto the ground.

Methods for cleaning up:

Absorb spillages with suitable absorbent material. Wash the area after the spillage is completely removed. Notification procedures:

Inform authorities if large amounts are involved.

Additional information:

N.A.

7. HANDLING AND STORAGE

Handling Local exhaust is recommended. Wear appropriate personal protection measures (see Section 8) to prevent skin or eye contact. Storage Store in closed original container in a dry place at -20°C. Keep away from food,

orage Store in closed original container in a dry place at -20°C. Keep away from food, drink and animal feeding stuffs.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with limit values that require monitoring at the workplace:

The exposure limit of DMSO (vapor and mist)is 250 ppm according to AIHA WEEL TWA.

General protective and hygienic measures:

The usual precautionary measures are the same as for chemicals handling.

Engineering measures:

Provide adequate ventilation.

Observe Occupational Exposure Limits and minimize the risk of inhalation of vapors. Personal protection:

- Respiratory protection: OSHA Standard 63FR 1152, January 8, 1998-approved respirator.
- Hand protection: Wear protective gloves. Be aware that the liquid may penetrate gloves.
- Frequent change is advisable. Suitable gloves can be recommended by the glove supplier. - Eye protection: Wear approved safety goggles in case risk of splashes.
- Eye protection. Wear approved safety goggles in case risk of splashes.
- Skin protection: Wear suitable protective clothing. Apron and long sleeves are recommended.



9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Liquid. Color: Orange Odor: Odorless or weak pH-value at 25°C: No data available. Boiling point: No data available. Melting point: No data available. Freezing point: No data available. Flash Point: No data available. Flammability / auto flammability: No data available. Explosive properties: No data available. Oxidizing properties: No data available. Vapor pressure: No data available. Relative density at 20°C: No data available. Solubility in/miscibility with water: Soluble Fat solubility: No data available. Vapor density: No data available. Evaporation rate: No data available. Conductivity: No data available. Viscosity: No data available.

10. STABILITY AND REACTIVITY

Stability: Stable under normal temperature conditions. Conditions to avoid: No conditions which may cause dangerous reactions. Materials to avoid: Strong oxidizing agents, strong bases. Hazardous decomposition products: At elevated temperature: Carbon dioxide, carbon monoxide, formaldehyde, sulfur oxides.

11. TOXICOLOGICAL INFORMATION

Chemical Name LD50 (dermal) LC50 (inhalation) LD50 (oral) No data available No data available No data available Novel Green Dye Dimethyl sulfoxide (DMSO) 14500 mg/kg (Rat) 40 g/kg (Rat) No data available Product information given on this section applies to concentrated solutions of Dimethyl sulfoxide or the substance in its pure form. Acute Effects: May cause skin and eye irritation. Chronic Effects: May cause injury to the liver and kidneys. Specific substance(s): Listed carcinogen: None.

12. ECOLOGICAL INFORMATION

*Information given on this section applies to concentrated solutions of Dimethyl sulfoxide or the substance in its pure form.

Ecotoxicity:

The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can make harmful or damaging effects on the environment.



Mobility: The product is water soluble and may spread in water systems. Persistence and degradability: No data available. Bioaccumulative potential: No data available on bioaccumulation. Other adverse effects: Not applicable.

13. DISPOSAL CONSIDERATIONS

Observe all national, state, and local regulations regarding disposal. Contact a licensed professional waste disposal service to dispose of this material.

14. TRANSPORT INFORMATION

This product does not contain hazardous compounds in quantities that require special handling precautions. Special precautions: None known. Proper shipping name: Not dangerous goods according to transport regulations. Transport regulations: No additional transport regulations are known to the supplier. IMDG: Non-hazardous for sea freight. ADR: Non-hazardous for road transport. RID: Non-hazardous for rail transport. ICAO/IATA: Non-hazardous for air transport.

15. REGULATORY INFORMATION

US

Labelling according to US Regulations: Caution: Avoid contact and inhalation.

US Regulatory information:

All intentional ingredients are listed on the TSCA Inventory.

All intentional ingredients are listed on the European EINECS Inventory.

16. OTHER INFORMATION

Hazard Ratings Health Hazard Fire Hazard Instability Special Hazard NFPA Λ n Hazard rating: 0 - Minimal; 1 - Slight; 2 - Moderate; 3 - Serious; 4 - Severe NFPA Label colored diamond code: Blue - Health; Red - Flammability; Yellow - Instability; White - Special Hazards Physical Hazard Health Hazard Flammability Personal Protection HMIS Hazard rating: 0 - Minimal; 1 - Slight; 2 - Moderate; 3 - Serious; 4 - Severe *- Chronic Health Effect

Personal Protection codes: B - Safety Glasses, Gloves

HMIS Label colored bar code: Blue - Health; Red - Flammability; Orange - Physical Hazards; White - Special The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The company shall not be held liable for any damage resulting from handling or from contact with the above product.