

Tray Plate

(#31001)

1. Description

The rectangle-shaped feature of the tray dishes provide larger surface area to ensure easy handling compared to round type dishes, for microbiological application.
The external dimensions of Tray plates are identical to standard multiwell plate for broader application.

2. Features

- ❖ Ideal for colony formation

3. Product Range

Cat.No.	Material	External Dimension w x l x h (mm)	Well Dimension w x l x h (mm)	Growth Area (cm ²)	Sterile	Packaging
31001	PS	127.94 x 85.5 x 16.25	109.64 x 73.46 x11.6	80.54	+	10 / 100

4. Certification


SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

7. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Plate		
Length: 85.5 mm Width: 127.94 mm Height: 9.00 mm	Length: 85.5 mm Width: 127.94 mm Height: 16.25 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 100

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-TP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 4, November 22th, 2018

