

Create Date: 2017-03-09 00:00:00/ Revision Date: 2017-08-01 00:00:00/ Print Date: 2019-10-29

SECTION 1: COMPANY AND PRODUCT INFORMATION

1.1 Product identifiers

Product name : PEG 6000 (Polyethylene glycol)
Product code : 41600211
CAS number : 25322-68-3
Synonyms : PEG

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, manufacture of substances

1.3 Details of the supplier of the safety data sheet

Address : Genelix International Inc, dba bioWORLD
4150 Tuller Rd. Suite 228
Dublin, OH 43017

Email : tech@bio-world.com

Phone : 614-792-8680, Toll free: 1-888-bio-PLUS

Fax : 614-792-8685

1.4 Emergency telephone number

Emergency phone : 1-888-bio-PLUS

SECTION 2: HAZARDS IDENTIFICATION

2.1 Classification of substance or mixture

Substance/mixture is not classified as hazardous.

2.2 GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Not classified as hazardous under the GHS.

2.3 Label elements and precautionary statements

Pictogram : Not available

Signal word : Not available

Hazard statement(s) : Not available

Precautionary statement(s) : No applicable precautionary statements.

2.4 Hazards not otherwise classified (HNOC) or not covered by GHS

No unclassified hazards known.

2.5 NFPA Rating

Health hazard : 0

Fire hazard : 0

Reactivity hazard : 0

2.6 HMIS Rating

Health hazard : 0
Chronic health hazard : 0
Reactivity hazard : 0
Flammability : 0
Physical hazard : 0

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substances

Substance	CAS#	EC#	Concentration
Polyethylene glycol M.F: (C ₂ H ₄ O) _n H ₂ O M.W: ~ 2000 g/mol	25322-68-3	500-038-2	95-100

3.2 Hazardous components & classification

Not available

SECTION 4: FIRST AID MEASURES

4.1 Description of first aid measures

General advice

Not available

In case of inhalation

If inhaled, move person to fresh air and monitor breathing. If not breathing, give artificial ventilation. Consult a physician if symptoms are severe or persistent.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water or eye wash solution as a precaution for 15 minutes. Consult a physician if symptoms are severe or persistent.

In case of ingestion

Never give anything by mouth to an unconscious person. Rinse mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

All known important symptoms are described in Section 2 and/or Section 11. No other important symptoms to report.

4.3 Indication of any immediate medical attention and special treatment needed

No special treatment indicated. Provide treatment in accordance with exhibited systems.

SECTION 5: FIREFIGHTING MEASURES

5.1 Suitable extinguishing media

Water spray, alcohol-resistant foam, dry chemical, and carbon dioxide are suitable

5.2 Unsuitable extinguishing media

Not available

5.3 Special hazards arising from the substance

Not available

5.4 Advice for firefighters

Wear protective gear, such as self-contained breathing apparatus, if necessary

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Avoid breathing vapors, mist or gas.

For personal protection see section 8

6.2 Environmental precautions

No special environmental precautions required.

6.3 Methods and material for containment and cleaning up

Sweep up or shovel. Keep in suitable, closed containers for disposal.

For proper disposal see section 13

SECTION 7: HANDLING AND STORAGE

7.1 Precautions for safe handling

Handle in accordance with good industrial hygiene and safety practice.

For precautions see section 2

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions : Store in cool, well ventilated area. Keep container tightly closed. Light Sensitive. Product is hygroscopic. Take precautions to avoid contact with atmospheric moisture. Store under Argon. Store at -20 °C

Incompatible materials : Strong oxidizing agents

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

8.2 Engineering controls

Showers, Eyewash stations, Ventilation systems

8.3 Personal protective equipment

Eye/face protection : Safety glasses with side-shields.

Skin protection : Wear protective gloves/clothing.

Body protection : Wear protective gloves/clothing.

Respiratory protection : If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Control of environmental exposure : No special environmental precautions required.

General hygiene considerations : Handle in accordance with good industrial hygiene and safety practice.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

- a) Appearance : White
- b) Physical states : Solid
- c) Odor : Not available
- d) Odor threshold : Not available
- e) Melting point : 4 - 8 °C (39 - 46 °F)
- f) Boiling point range : Not available
- g) pH : Not available
- h) Density : Not available
- i) Flash point : Not available
- j) Evaporation rate : Not available
- k) Flammability : Not available
- l) Upper/lower flammability or explosive limits: : Not available
- m) Vapor pressure : < 0.01 hPa (< 0.01 mmHg) at 20 °C (68 °F)
- n) Vapor density : Not available
- o) Relative density : 1.128 g/mL
- p) Water solubility : Not available
- q) Partition coefficient:n-octanol/water : Not available
- r) Autoignition temperature : Not available
- s) Decomposition temperature : > 300 °C (> 572 °F) -
- t) Kinematic viscosity : Not available
- u) Explosive properties : Not available
- v) Oxidizing properties : Not available
- w) Solubility in other solvents : Not available
- x) Surface tension : Not available

SECTION 10: STABILITY AND REACTIVITY

10.1 Reactivity

No special reactivity is known. Non-reactive under normal use.

10.2 Chemical stability

Product is stable when stored as recommended.

10.3 Stability note(s)

No special or unusual instability known.

10.4 Polymerization

No known polymerization possible.

10.5 Possibility of hazardous reactions

No known hazardous reactions are possible.

10.6 Incompatible materials

Strong oxidizing agents.

10.7 Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1 Acute toxicity

LD50 Oral : No toxicity data available.

LD50 Dermal : No toxicity data available.

LC50 Inhalation : No toxicity data available.

11.2 Skin corrosion/irritation

No skin/corrosion irritation data available.

11.3 Serious eye damage/eye irritation

No eye damage/irritation data available.

11.4 Respiratory or skin sensitization

No sensitization data available.

11.5 Germ cell mutagenicity

No mutagenicity data available.

11.6 Carcinogenicity

IARC : No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH : No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP : No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA : No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

11.7 Reproductive toxicity

No reproductive toxicity data available.

11.8 Specific target organ toxicity – single exposure

No specific organ toxicity data available.

11.9 Specific target organ toxicity – repeated exposure

No specific organ toxicity data available.

11.10 Aspiration hazard

No aspiration hazard data available.

11.11 Additional Information

RTECS: No RTECS data available.

SECTION 12: ECOLOGICAL INFORMATION

12.1 Toxicity

No ecological toxicity data available.

12.2 Aquatic toxicity

No aquatic toxicity data available.

12.3 Persistence and degradability

No persistence/degradability data available.

12.4 Bioaccumulative potential

No bioaccumulation data available.

12.5 Mobility in soil

No soil mobility data available.

12.6 Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment is not required/was not conducted.

12.7 Other adverse effect

No other adverse effect data available.

SECTION 13: DISPOSAL CONSIDERATIONS

13.1 Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

13.2 Packaging

Dispose of as unused product.

13.3 Recommendation

Disposal must be made according to official regulations.

SECTION 14: TRANSPORTATION INFORMATION

14.1 DOT (US)

Not dangerous goods

14.2 IMDG

Not dangerous goods

14.3 IATA

Not dangerous goods

SECTION 15: REGULATORY INFORMATION

15.1 SARA

SARA 302 Components

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

15.2 Clean water act (CWA)

No chemicals are present in this product that are subject to regulation under the Clean Water Act.

15.3 Right to know components

Massachusetts : No chemicals are present which require disclosure under the Massachusetts Right to Know Act.

Pennsylvania : 5-Hydro-9-hydroxy-poly(oxy-1,2-ethanediyl), n ~400 CAS-No-25322-68-3

New Jersey : 5-Hydro-9-hydroxy-poly(oxy-1,2-ethanediyl), n ~400 CAS-No-25322-68-3

California : This product contains no chemicals which are known to the State of
proposition 65 California to cause cancer, or birth defects or other reproductive harm.
components

SECTION 16: OTHER INFORMATION

16.1 Disclaimer

This product is offered by bioWORLD for research, laboratory or further manufacturing use. Not for human use or consumption. The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchant-ability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall bio-WORLD.com be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if bio-WORLD.com has been advised of the possibility of such damages.

16.2 Copyright

Copyright 2019 www.bio-WORLD.com. License granted to make paper copies for internal use only.

16.2 Preparation Information

Create Date: 2017-03-09 00:00:00/ Revision Date: 2017-08-01 00:00:00/ Print Date: 2019-10-29