

Square Dish (#10125)

1. Description

Square Dishes have been widely used for microbiology application such as colony tracing and picking. Large surface area with gridded bottoms of the Square Dishes are useful for tracing the location of individual colony.

2. Features

- ❖ Ideal for colony count or determination of plating efficiency
- ❖ Useful for counting or colony picking

3. Product Range

Cat.No.	Material	External Dimension w x l x h (mm)	Well Dimension w x l x h (mm)	Growth Area (cm ²)	Sterile	Packaging
11125	PS	126.40 x 126.40 x 20.00	118.70 x 118.70 x 11.60	139.00	+	10 / 240
10245	PS	243.00 x 243.00 x 19.30	227.50 x 227.50 x 13.20	500.00	+	5 / 20
11245	PS	243.00 x 243.00 x 27.30	227.50 x 227.50 x 21.20	500.00	+	5 / 20

4. Certification


SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Specifications & Technical Drawing

(1) Square Dish (# 11125)
- Specifications

Dimensions		Material (Resin)	Packaging
Length:	126.40 mm	Crystal-grade polystyrene (PS)	▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 240
Width:	126.40 mm		
Height:	20.00 mm		


- Technical Drawing



(2) Square Dish (# 10245)
- Specifications

Dimensions		Material (Resin)	Packaging
Length:	243.00 mm	Crystal-grade polystyrene (PS)	▪ Pieces/Sleeve: 5 ▪ Pieces/Case: 20
Width:	243.00 mm		
Height:	19.30 mm		


- Technical Drawing



(3) Square Dish (# 11245)
 - Specifications

Dimensions		Material (Resin)	Packaging
Length:	243.00 mm	Crystal-grade polystyrene (PS)	• Pieces/Sleeve: 5 • Pieces/Case: 20
Width:	243.00 mm		
Height:	27.30 mm		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr