

250 ml Conical Tube (#558250)

1. Description

250 ml Conical tube can be used in storage of samples and centrifugation of large volume samples, between 50 ml to 250 ml. The product is provided with a dedicated Styrofoam rack for user convenience.

2. Features

- ❖ Tight sealing screw cap
- ❖ External graduations with marking area
- ❖ Packed in sleeves
- ❖ Polypropylene conical tubes
- ❖ Provided with styrofoam rack

3. Functional Test

Non-cytotoxic	:	Passed
Non-pyrogenic	:	Passed
DNase / RNase – free	:	Passed
Human DNA – free	:	Passed

4. Product Range

In Sleeves

Cat.No.	Bottom Type	Material (Body/Cap)	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
558250	Conical	PP / PP	60.00 x 148.00	250.00	5,000 x g	+	6 in sleeve / 54

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs




6. Specifications & Technical Drawing

- Specifications

Dimensions			Material (Resin)	Packaging
Cap	Diameter:	42.50 mm	Polypropylene (PP)	<div>▪ Pieces/sleeve: 6</div> <div>▪ Pieces/case: 54</div>
	Height:	15.00 mm		
Body	Diameter:	60.00 mm	Polypropylene(PP)	
	Height:	147.00 mm		
	Volume:	250.00 ml		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr