

25ml Snap Tube (#50425/54425)

1. Description

25 ml Snap tubes conical type are widely used in laboratories for multiple purposes, from storage of samples to high speed centrifugation. Snap tubes can assist users during experiments by enabling easy one-hand opening and closing of the tubes by providing snap type cap.

2. Features

- ❖ Snap cap with tight sealing
- ❖ Imprinted external graduations with marking area
- ❖ Packed in zipper bags
- ❖ 25 ml Snap Tube Adapter for 50 ml rotor (Cat.No.52025)

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase Free	:	Passed
DNA free	:	Passed

4. Product Range

Cat.No.	Material	Color	Bottom Type	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
50425	PP	Clear	Conical	28.30 x 77.15	25.00	18,000 x g	+	25 in sleeve / 200
54425	PP	Amber	Conical	28.30 x 77.15	25.00	18,000 x g	+	25 in sleeve / 200

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.


6. Specifications & Technical Drawing

(1) 25 ml Snap Tube (#50425/54425)

- Specifications

Dimensions		Material (Resin)	Packaging
Cap	Diameter: 34.10 mm Height: 7.50 mm	Polypropylene(PP)	<ul style="list-style-type: none"> ▪ Pieces/sleeve:25 ▪ Pieces/case: 200
Body	Diameter: 28.54 mm Height: 77.15 mm Volume: 25.00 ml	Polypropylene(PP)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr