

SPLAmber™ 5mL Snap Tube

(#54005/54105)

1. Description

SPLAmber™ 5mL Snap Tubes are suitable for storage of light-sensitive samples and high speed centrifugation. Materials of SPLAmber™ significantly decrease light transmission over the visible wavelength spectrum compared to clear polypropylene. Snap type cap can assist users during experiments by enabling easy one-hand opening and closing of the tubes..

2. Features

- ❖ Light protection for light-sensitive samples
- ❖ Snap cap with tight sealing
- ❖ Sterile (Cat. No.54105)
- ❖ 5ml Snap Tube Adapter for 15ml rotor (Cat.No.52005)

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase Free	:	Passed
DNA free	:	Passed

4. Product Range

Cat.No.	Material	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
54005	PP	20.00 x 62.00	5	25,000 x g	-	100 in sleeve / 200
54105	PP	20.00 x 62.00	5	25,000 x g	+	100 in sleeve / 200

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs




SPL-TDS-A5ST

7. Specifications & Technical Drawing

- Specifications

Dimensions			Material (Resin)	Packaging
Cap	Diameter: Height:	22.08 mm 7.00 mm	Polypropylene(PP)	▪ Pieces/sleeve:100 ▪ Pieces/case: 200
Body	Diameter: Height: Volume:	15.80 mm 59.30 mm 5.00 ml	Polypropylene(PP)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-A5ST

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, July 19th, 2016

