

## Agarose Gel Scalpel (# 90080)

### 1. Description

Agarose Gel Scalpel can be easily used to extract DNA or RNA bands from agarose gel after electrophoresis. The plastic blade makes it safer than a typical knife and allows you to conduct experiments without exposing your hand to transilluminator screens.

### 2. Features

- ❖ Structure capable of cutting gel avoiding UV transilluminator screens
- ❖ Sharp end: can be used to transfer excised piece of agarose gel
- ❖ Easy cutting
- ❖ Disposable

### 3. Functional Test

- ❖ RNase/DNase Free
- ❖ Human-DNA Free

### 4. Product Range

| Cat.No. | Materials | Total Length (mm) | Sterile | Packaging |
|---------|-----------|-------------------|---------|-----------|
| 90080   | PS        | 153.54            | -       | 10 / 20   |

### 5. Product Photographs



### 6. Certification


SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

7. Specifications & Technical Drawing

(1) Agarose Gel Scalpel (# 90080)  
 - Specifications

| Body Dimensions |           | Materials | Sterile | Packaging              |
|-----------------|-----------|-----------|---------|------------------------|
| Length:         | 153.54 mm | PS        | -       | ▪ Pieces/Case: 10 / 20 |
| Width:          | 13.84 mm  |           |         |                        |
| Height:         | 4.60 mm   |           |         |                        |

- Technical Drawing



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:  
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)  
 To place an order, contact your local distributor or  
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-Gel Extractor

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
 © SPL Life Sciences Co., Ltd., Version 1, SEP. 4<sup>nd</sup>, 2020

